


MONDRIAN™

TRANSFORAMINAL LUMBAR
INTERBODY FUSION DEVICES

MONDRIAN™ transforaminal lumbar cages are titanium plated PEEK lumbar interbody fusion implants designed to accommodate a wide range of patient anatomies and surgical preferences.




Large graft window allows for maximum graft volume and greater likelihood of fusion

Curved design mimics patients' vertebral anatomy

Aggressive serrations ensure endplate fixation

Self distracting tip for ease of insertion

Lateral apertures for improved visualization via imaging



CTL Amedica's proprietary **TiCro®** design offers significantly greater surface area, improving endplate contact.

This unique surface geometry enhances bone interlocking properties and helps to ensure cage placement.



CTL AMEDICA
rethink what's possible®

MONDRIAN™ TLIF Curved Ti-PEEK Cage: W11 mm x L28mm

	0 Degree*	9 Degree	15 Degree
Height (mm)	Part Number	Part Number	Part Number
7	117.2637		
8	117.2638	117.2718	
9	117.2639	117.2719	117.2799
10	117.2640	117.2720	117.2800
11	117.2641	117.2721	117.2801
12	117.2642	117.2722	117.2802
13	117.2643	117.2723	117.2803
14	117.2644	117.2724	117.2804
15	117.2645	117.2725	117.2805

*available upon request

MONDRIAN™ TLIF Curved Ti-PEEK Cage: W11 mm x L32mm

	0 Degree*	9 Degree	15 Degree
Height (mm)	Part Number	Part Number	Part Number
7	117.2657		
8	117.2658	117.2738	
9	117.2659	117.2739	117.2819
10	117.2660	117.2740	117.2820
11	117.2661	117.2741	117.2821
12	117.2662	117.2742	117.2822
13	117.2663	117.2743	117.2823
14	117.2664	117.2744	117.2824
15	117.2665	117.2745	117.2825

*available upon request



CTL AMERICA
rethink what's possible®